DE: robomove® expander vh 2,0 m
cos101659va02

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EN: robomove® expander vh 2.0 m
cos101659va02

Elastic cable system for vertical or horizontal wall mount

Included in delivery:

2 nooses for hand or leg

8 screws and dowels for wall mount

Instructions for use and warnings:

Before using the robomove® expander it has to be mounted on a wall. Considering a safety factor 5 the mounting has to bear a pulling load of max. 100 kg. So the mounting of the robomove® expander has to be performed by qualified personnel only. The added mounting material is not applicable on every wall. So the decision about the right mounting material has to be performed by a specialist also.

The h/p/cosmos robomove® expander is an elastic cable system for wall mount and therewith compatible to all kind of treadmill models.

The patient's arms/legs are attached to the elastic cables with the nooses. The elastic cables generate movement support or resistance.

By adjusting the angle of support/resistance movement correction is possible.

The adjustment is very fast and simple.

Scales give reproducibility of the adjustments.

The patented tension adjustment module consists of continuously variable rubber cable pulls with scales for load display.

The traction device is based on the patented Bodyspider Technology. See patents EP1221331 and WO9823334.

Examples of use include:

Motion support, mobilization, locomotion, gait training and gait correction for orthopaedic or neurological problems, coordination and functional training, strength and endurance training.

material:

stainless steel tube 30x30x2

measurement device:

100x200x2000 mm

measurement box:

300x300x2100 mm

weight device:

7.0 kg

weight incl. box:

8.0 kg

surface slide:

powder coating ral9007 grey aluminium

surface tube:

electro polished

expander length:

6.0 m

max. pulling load/force expanders:

grey 5 kg (50 N), black 4 kg (40 N) approx.

Warnings!

When using the robomove® expander in combination with a treadmill, only an h/p/cosmos treadmill model may be used. Furthermore a safety system to prevent from falling is obligatory when using in combination with an h/p/cosmos treadmill. We recommend the h/p/cosmos safety arch with chest belt, rope and automatic stop of the treadmill or an unweighting system that prevents from falling. The exercise has to be permanently supervised by a physiotherapist or a medical doctor in the field of orthopaedics or neurology.

Disregard of warnings and/or the use of the robomove® expander on a treadmill without a fall protection system is strictly prohibited and may lead to accidents and injuries with fatal consequences.

All kinds of installation, repair work and maintenance are to be only performed by trained and authorized technicians who have been certified by h/p/cosmos.

Various patents pending.

robomove® is a registered and protected trade mark of Franz Harrer.

Made in Germany for h/p/cosmos sports & medical gmbh

Empfohlene Montageposition / suggested mounting position: