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Note:

- Determine the length of the running belt when installed and correctly tensioned. **Running belt length in mm:** _____
- Switch on the treadmill and connect it to the software.
[h/p/cosmos para control® 4.2 \(up to MCU5\)](#)
[h/p/cosmos para control® 6 \(from MCU6\)](#)
- The speed unit must be set to km/h.
- Set the target speed via the software (see table).
- When the target speed is reached (display value is constant and does not flash), determine the time required for exactly 10 belt revolutions and enter it here:

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Calculation (example):

Speed = running belt length / revolutions / time x 0.36

8.09 km/h = 3348 / 10 / 14.90 x 0.36

The columns "2. Check" and "3. Check" are comparison values and must be checked in any case.

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Note:

The speed control is performed without rotor and no load. For this reason, it is possible that the calculated values are slightly higher for devices without speed measurement than the indicated value. Under load, the value falls then back slightly, and thus corresponds to the display.

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* If the distance X at 0% elevation will not be the same as mentioned, the levelling sockets must be adjusted according to this distance. Otherwise, the machine is not levelled correctly. For machines with frame type B, the values on the left and right side should be similar, to ensure the synchronization of the elevation system. The differences between the left and right side must be smaller than 6mm. If more, please contact h/p/cosmos service department immediately. If the measured value is not within the limits, a troubleshooting to find the source is necessary. All abnormalities must be noted on the service report.

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