



America

# CERTIFICATE

No. U8 13 02 45283 009

**Holder of Certificate:** h/p/cosmos  
sports & medical gmbhAm Sportplatz 8  
83365 Nussdorf - Traunstein  
GERMANY**Production Facility(ies):** 45283**Certification Mark:****Product:** Therapeutic training devices  
Treadmill Ergometer (running machine)**Model(s):** 1) Product Family: treadmill h/p/ cosmos 150/50  
2) Product Family: treadmill h/p/ cosmos 170-190/65  
for more information please see attachment**Parameters:** Protection Class: I  
Applied Part: B  
  
for more information please see attachment**Remarks:**  
The certificate is valid for USA and Canada.  
Other countries might have other requirements,  
which were not part of this certification.**Tested according to:** UL 60601-1:2003  
CAN/CSA C22.2 No.601.1-M90

The product was voluntarily tested according to the relevant safety requirements noted above. It can be marked with the certification mark above. The mark must not be altered in anyway. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC Guide 67. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited certification body.

**Test report no.:** 713003265A**Date,** 2013-03-01

Page 1 of 3





Attachment to NRTL Certificate no. U8 13 02 45283 009

Factories: 45823

Product: Therapeutic training devices  
Treadmill Ergometer (running machine)

Model(s): 1) Product family: treadmill h/p/cosmos 150 / 50 consisting of the following models:

article no.	variant	variant model name
cos30000	va05	h/p/cosmos stratos lt med
cos30000	va06	h/p/cosmos stratos med
cos30000	va07	h/p/cosmos mercury lt med
cos30000	va08	h/p/cosmos mercury med
cos30001	va01	h/p/cosmos locomotion 150/50 E med
cos30001-01	va02	h/p/cosmos locomotion 150/50 DE med
cos30002	va01	h/p/cosmos Gaitway II F
cos30002	va02	h/p/cosmos Gaitway II S
cos30000	va09	CareFusion LE 100CE 150/50
cos30000	va10	CareFusion LE 200CE 150/50
cos30000	va15	COSMED T150 MED
cos30000	va16	COSMED T150 D MED
cos30000	va17	COSMED T150 E MED
cos30000	va18	COSMED T150 DE MED
cos30000	va19	proxomed Kardiomed Mill F
cos30000	va20	proxomed Kardiomed Mill S
cos30000	va23	SCHILLER MTM 1500 med
cos30000	va28	LMT ProMedPlus 150/50

2) Product family: treadmill h/p/cosmos 170-190 / 65 consisting of the following models:

article no.	variant	variant model name
cos30003	va17	h/p/cosmos stellar lt med
cos30003	va18	h/p/cosmos stellar med
cos30003	va19	h/p/cosmos quasar lt med
cos30003	va20	h/p/cosmos quasar med
cos30003	va01	CareFusion LE 250CE 170/65
cos30003	va02	CareFusion LE 300CE 170/65
cos30003	va03	CareFusion LE 500CE 170/65



cos30003	va08	COSMED T170 MED
cos30003	va09	COSMED T170 D MED
cos30003	va10	COSMED T170 E MED
cos30003	va11	COSMED T170 DE MED
cos30003	va12	COSMED T170 DE sportmed
cos30004	va01	h/p/cosmos pulsar lt
cos30004	va03	h/p/cosmos pulsar
cos30024	va01	h/p/cosmos locomotion 190/65 E med
cos30024	va03	h/p/cosmos locomotion 190/65 DE med

Parameters: Rated Input:  
 1) 100 / 110 / 200 / 208 / 230 V AC  
 2) 100 / 110 / 200 / 208 / 230 V AC

Rated Frequency:  
 1) 50 / 60 Hz  
 2) 50 / 60 Hz

Rated Input Power long time:  
 1) 1500 VA  
 2) 1600 VA

Rated Input Power momentary:  
 1) 3400 VA  
 2) 3500 VA

Remarks: The certificate is valid for USA and Canada.  
 Other countries might have other requirements,  
 Which were not part of this certification.

Date, 2013-03-01