



torqualizer® 600ef med



torqualizer 600ef med

In addition to the treadmill, the bicycle ergometer is a classic machine in cardio training and another ideal test- and training device for the professional user. Independent from power supply (no power cable!) the unit can be assured even with most of the various options, which guarantees flexibility and saves money. The optional length adjustable crank arm is an extremely important feature for all users with limitations in the range of motion.

Functions & Features

- easy to handle, multicolor, backlit display with Bio-Feedback
- easy step through
- horizontally and vertically adjustable saddle
- coded Polar® pulse receiver
- RS232 interface

Options & order number



ergometer torqualizer 600ef med	cos30021ef-med600
power upgrade 750 watts	cos30021ef-pup
therapy program package: 8 more programs	cos30021ef-thp
test program package	cos30021ef-tep
length adjustable crank arms	cos30021ef-acr
RFID-System / wireless package	cos30021ef-rfi
ear pulse system	cos30021ef-eps
increase of max. user weight to 200 kg	cos30021ef-muw
packing in crate	cos30021ef-pic



Specifications

- grid independent, no power cord
- calibrated drive / brake unit
- hybrid-brake-system (HBS)
- mono belt drive, self-adjusting and silent
- power range: 15-500 watts (optionally 750 watts) rpm independent, 15 - 1000 watts rpm dependent workload steps: 5 watts
- speed range: 20 – 120 rpm/min
- manufactured according to: DIN EN ISO 20957-1 SA, DIN EN ISO 20957-5 SA, DIN EN 60601-1, DIN EN 60601-1-2
- CE based on the guidelines for medical devices according to directive 93/42 EEC, class IIa
- max user weight: 150 kg
- Size (LxWxH): 120 x 65 x 155 cm
- Weight: 63 kg unpacked
- Made in Germany

dealer:

contact

h/p/cosmos sports & medical gmbh
 Am Sportplatz 8
 83365 Nussdorf-Traunstein
 Germany

phone: +49 86 69 86 42 0
 fax: +49 86 69 86 42 49

sales@h-p-cosmos.com
 www.h-p-cosmos.com

skype: @h-p-cosmos.com (search & select name)
 youtube: youtube.com/hpcosmos
 twitter: twitter.com/hpcosmos
 facebook: facebook.com/hpcosmos