

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT
(IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Name and address of the applicant

Name and address of the manufacturer

Name and address of the factory

Ratings and principal characteristics

Trade mark (if any)

Model/type Ref.

Additional information (if necessary)

A sample of the product was tested and found to be in conformity with

as shown in the Test Report Ref. No. which forms part of this certificate

Treadmill-Ergometer

h/p/cosmos

sports & medical gmbh

Am Sportplatz 8

83365 Nussdorf - Traunstein

GERMANY

h/p/cosmos sports & medical gmbh

Am Sportplatz 8, 83365 Nussdorf - Traunstein, GERMANY

h/p/cosmos sports & medical gmbh

Am Sportplatz 8, 83365 Nussdorf - Traunstein, GERMANY

Rated voltage: 200 V – 240 V_{AC} or 120 V_{AC}
 Rated frequency: 50 – 60 Hz
 Rated current: long time 6,0 A & momentary 13,5 A
 or
 long time 9.5 A & momentary 23.5 A
 Protection class: I
 Type of applied part: B
 Software version: MCU 5, v1.08.02

h/p/cosmos

pluto[®] med, pluto[®] It med, mercury med, motion sprint 600 PL med

Certificate DE 3 - 40682 issued 2017-01-20 is replaced by this version due to technical changes

IEC 60601-1:2005
 IEC 60601-1:2005/AMD1:2012
 IEC 60601-1-6:2010
 IEC 60601-1-6:2010/AMD1:2013

713141002

This CB Test Certificate is issued by the National Certification Body

CB 045283 0024 Rev. 00

Date, 2019-07-02



(Paul Bielawski)



Page 1 of 1

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany

Product Service