

model name product code / article number	name plate / type plate with product family name, model name and technical data	name plate / type plate with UDI PI, manufacturer ID, serial number and manufacturing date (date is example only and will be replaced by actual date)																																								
pluto® med cos30026-01va02	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td colspan="3">pluto® med</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.1</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 250 kg / 551 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00	model:	pluto® med			U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1			current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 250 kg / 551 lbs				<table><tr><td>(21)cos30026-01va02-0001 (11)181218</td><td></td><td>2018-12-18</td><td></td></tr><tr><td></td><td></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI (01)4050588002336</td><td></td><td>SN cos30026-01va02-0001</td><td></td></tr><tr><td colspan="4">Made in Germany</td></tr></table>	(21)cos30026-01va02-0001 (11)181218		2018-12-18				h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI (01)4050588002336		SN cos30026-01va02-0001		Made in Germany			
product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00																																					
model:	pluto® med																																									
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1																																							
current input: long time 6.0 A / momentary 13.5 A																																										
max weight on running surface: 250 kg / 551 lbs																																										
(21)cos30026-01va02-0001 (11)181218		2018-12-18																																								
		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																								
UDI (01)4050588002336		SN cos30026-01va02-0001																																								
Made in Germany																																										
pluto® med OEM cos30026-01va04	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td colspan="3">pluto® med OEM</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.1</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 250 kg / 551 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00	model:	pluto® med OEM			U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1			current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 250 kg / 551 lbs				<table><tr><td>(21)cos30026-01va14-0001 (11)210614</td><td></td><td>2021-06-14</td><td></td></tr><tr><td></td><td></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI (01)4050588003173</td><td></td><td>SN cos30026-01va14-0001</td><td></td></tr><tr><td colspan="4">Made in Germany</td></tr></table>	(21)cos30026-01va14-0001 (11)210614		2021-06-14				h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI (01)4050588003173		SN cos30026-01va14-0001		Made in Germany			
product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00																																					
model:	pluto® med OEM																																									
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1																																							
current input: long time 6.0 A / momentary 13.5 A																																										
max weight on running surface: 250 kg / 551 lbs																																										
(21)cos30026-01va14-0001 (11)210614		2021-06-14																																								
		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																								
UDI (01)4050588003173		SN cos30026-01va14-0001																																								
Made in Germany																																										
pluto® lt med cos30027-01va02	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td colspan="3">pluto® lt med</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.1</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 250 kg / 551 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00	model:	pluto® lt med			U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1			current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 250 kg / 551 lbs				<table><tr><td>(21)cos30027-01va02-0001 (11)180817</td><td></td><td>2018-08-17</td><td></td></tr><tr><td></td><td></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI (01)4050588002435</td><td></td><td>SN cos30027-01va02-0001</td><td></td></tr><tr><td colspan="4">Made in Germany</td></tr></table>	(21)cos30027-01va02-0001 (11)180817		2018-08-17				h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI (01)4050588002435		SN cos30027-01va02-0001		Made in Germany			
product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00																																					
model:	pluto® lt med																																									
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1																																							
current input: long time 6.0 A / momentary 13.5 A																																										
max weight on running surface: 250 kg / 551 lbs																																										
(21)cos30027-01va02-0001 (11)180817		2018-08-17																																								
		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																								
UDI (01)4050588002435		SN cos30027-01va02-0001																																								
Made in Germany																																										
pluto® ef med cos30027-01va14	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td colspan="3">pluto® ef med</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.1</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 150 kg / 330 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00	model:	pluto® ef med			U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1			current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 150 kg / 330 lbs				<table><tr><td>(21)cos30027-01va14-0001 (11)200403</td><td></td><td>2020-04-03</td><td></td></tr><tr><td></td><td></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI (01)4050588002565</td><td></td><td>SN cos30027-01va14-0001</td><td></td></tr><tr><td colspan="4">Made in Germany</td></tr></table>	(21)cos30027-01va14-0001 (11)200403		2020-04-03				h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI (01)4050588002565		SN cos30027-01va14-0001		Made in Germany			
product family:	treadmill h/p/cosmos 150/50			CE 0123 MD	IP00																																					
model:	pluto® ef med																																									
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1																																							
current input: long time 6.0 A / momentary 13.5 A																																										
max weight on running surface: 150 kg / 330 lbs																																										
(21)cos30027-01va14-0001 (11)200403		2020-04-03																																								
		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																								
UDI (01)4050588002565		SN cos30027-01va14-0001																																								
Made in Germany																																										

pluto® lt med OEM cos30027-01va16	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 </td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td>pluto® lt med OEM</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.1</td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 250 kg / 551 lbs</td></tr></table>  	product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00	model:	pluto® lt med OEM	U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1	current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 250 kg / 551 lbs				<table><tr><td>(21)cos30027-01va16-0001</td><td>2022-02-28</td><td>h/p/cosmos</td></tr><tr><td>(11)220228</td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI</td><td>SN</td><td>cos30027-01va16-0001</td></tr><tr><td>(01)4050588003210</td><td colspan="2">Made in Germany</td></tr></table> 	(21)cos30027-01va16-0001	2022-02-28	h/p/cosmos	(11)220228	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI	SN	cos30027-01va16-0001	(01)4050588003210	Made in Germany	
product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00																													
model:	pluto® lt med OEM																															
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.1																													
current input: long time 6.0 A / momentary 13.5 A																																
max weight on running surface: 250 kg / 551 lbs																																
(21)cos30027-01va16-0001	2022-02-28	h/p/cosmos																														
(11)220228	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
UDI	SN	cos30027-01va16-0001																														
(01)4050588003210	Made in Germany																															
mercury® med cos30000-02va02	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 </td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td>mercury® med</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.4</td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 661 lbs</td></tr></table>  	product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00	model:	mercury® med	U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4	current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 300 kg / 661 lbs				<table><tr><td>(21)cos30000-02va02-0001</td><td>2018-12-17</td><td>h/p/cosmos</td></tr><tr><td>(11)181217</td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI</td><td>SN</td><td>cos30000-02va02-0001</td></tr><tr><td>(01)4050588002602</td><td colspan="2">Made in Germany</td></tr></table> 	(21)cos30000-02va02-0001	2018-12-17	h/p/cosmos	(11)181217	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI	SN	cos30000-02va02-0001	(01)4050588002602	Made in Germany	
product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00																													
model:	mercury® med																															
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4																													
current input: long time 6.0 A / momentary 13.5 A																																
max weight on running surface: 300 kg / 661 lbs																																
(21)cos30000-02va02-0001	2018-12-17	h/p/cosmos																														
(11)181217	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
UDI	SN	cos30000-02va02-0001																														
(01)4050588002602	Made in Germany																															
stratos® med cos30000-02va10	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 </td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td>stratos® med</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.4</td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 661 lbs</td></tr></table>  	product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00	model:	stratos® med	U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4	current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 300 kg / 661 lbs				<table><tr><td>(21)cos30000-02va10-0001</td><td>2020-10-05</td><td>h/p/cosmos</td></tr><tr><td>(11)201005</td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI</td><td>SN</td><td>cos30000-02va10-0001</td></tr><tr><td>(01)4050588002893</td><td colspan="2">Made in Germany</td></tr></table> 	(21)cos30000-02va10-0001	2020-10-05	h/p/cosmos	(11)201005	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI	SN	cos30000-02va10-0001	(01)4050588002893	Made in Germany	
product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00																													
model:	stratos® med																															
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4																													
current input: long time 6.0 A / momentary 13.5 A																																
max weight on running surface: 300 kg / 661 lbs																																
(21)cos30000-02va10-0001	2020-10-05	h/p/cosmos																														
(11)201005	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
UDI	SN	cos30000-02va10-0001																														
(01)4050588002893	Made in Germany																															
mercury® med OEM cos30000-02va12	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 </td><td rowspan="2">IP00</td></tr><tr><td>model:</td><td>mercury® med OEM</td></tr><tr><td>U: 200 V - 240 V ~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 8.4</td></tr><tr><td colspan="4">current input: long time 6.0 A / momentary 13.5 A</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 661 lbs</td></tr></table>  	product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00	model:	mercury® med OEM	U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4	current input: long time 6.0 A / momentary 13.5 A				max weight on running surface: 300 kg / 661 lbs				<table><tr><td>(21)cos30000-02va12-0001</td><td>2021-06-14</td><td>h/p/cosmos</td></tr><tr><td>(11)210614</td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI</td><td>SN</td><td>cos30000-02va12-0001</td></tr><tr><td>(01)4050588003197</td><td colspan="2">Made in Germany</td></tr></table> 	(21)cos30000-02va12-0001	2021-06-14	h/p/cosmos	(11)210614	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI	SN	cos30000-02va12-0001	(01)4050588003197	Made in Germany	
product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP00																													
model:	mercury® med OEM																															
U: 200 V - 240 V ~	f: 50 Hz - 60 Hz	class: S, I, A	code: 8.4																													
current input: long time 6.0 A / momentary 13.5 A																																
max weight on running surface: 300 kg / 661 lbs																																
(21)cos30000-02va12-0001	2021-06-14	h/p/cosmos																														
(11)210614	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
UDI	SN	cos30000-02va12-0001																														
(01)4050588003197	Made in Germany																															
locomotion® 150/50 DE med cos30001-01va02	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 150/50</td><td rowspan="2">CE 0123 </td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>locomotion® 150/50 DE med</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td>code: 1.5</td></tr><tr><td colspan="4">power input: long time 1500 VA / momentary 3400 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>  	product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP20	model:	locomotion® 150/50 DE med	U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A	code: 1.5	power input: long time 1500 VA / momentary 3400 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30001-01va02-0001</td><td>2015-11-11</td><td>h/p/cosmos</td></tr><tr><td>(11)151111</td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>UDI</td><td>SN</td><td>cos30001-01va02-0001</td></tr><tr><td>(01)4050588001506</td><td colspan="2">Made in Germany</td></tr></table> 	(21)cos30001-01va02-0001	2015-11-11	h/p/cosmos	(11)151111	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		UDI	SN	cos30001-01va02-0001	(01)4050588001506	Made in Germany	
product family:	treadmill h/p/cosmos 150/50	CE 0123  	IP20																													
model:	locomotion® 150/50 DE med																															
U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A	code: 1.5																													
power input: long time 1500 VA / momentary 3400 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30001-01va02-0001	2015-11-11	h/p/cosmos																														
(11)151111	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
UDI	SN	cos30001-01va02-0001																														
(01)4050588001506	Made in Germany																															

stellar® lt med cos30003va17	<table><tr><td>product family:</td><td colspan="2">treadmill h/p/cosmos 170-190/65</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="2">stellar® lt med</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="3">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20	model:	stellar® lt med		U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 1600 VA / momentary 3500 VA			max weight on running surface: 300 kg / 660 lbs			<table><tr><td>(21)cos30003va17-0001 (11)151113</td><td> 2015-11-13</td><td></td></tr><tr><td> </td><td> h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td rowspan="2"></td></tr><tr><td>(01)4050588000493</td><td> cos30003va17-0001 Made in Germany</td></tr></table>	(21)cos30003va17-0001 (11)151113	 2015-11-13		 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588000493	 cos30003va17-0001 Made in Germany
product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20																									
model:	stellar® lt med																												
U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A																											
power input: long time 1600 VA / momentary 3500 VA																													
max weight on running surface: 300 kg / 660 lbs																													
(21)cos30003va17-0001 (11)151113	 2015-11-13																												
 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																												
(01)4050588000493	 cos30003va17-0001 Made in Germany																												
stellar® med cos30003va18	<table><tr><td>product family:</td><td colspan="2">treadmill h/p/cosmos 170-190/65</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="2">stellar® med</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="3">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20	model:	stellar® med		U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 1600 VA / momentary 3500 VA			max weight on running surface: 300 kg / 660 lbs			<table><tr><td>(21)cos30003va18-0001 (11)151113</td><td> 2015-11-13</td><td></td></tr><tr><td> </td><td> h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td rowspan="2"></td></tr><tr><td>(01)4050588000509</td><td> cos30003va18-0001 Made in Germany</td></tr></table>	(21)cos30003va18-0001 (11)151113	 2015-11-13		 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588000509	 cos30003va18-0001 Made in Germany
product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20																									
model:	stellar® med																												
U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A																											
power input: long time 1600 VA / momentary 3500 VA																													
max weight on running surface: 300 kg / 660 lbs																													
(21)cos30003va18-0001 (11)151113	 2015-11-13																												
 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																												
(01)4050588000509	 cos30003va18-0001 Made in Germany																												
quasar® lt med cos30003va19	<table><tr><td>product family:</td><td colspan="2">treadmill h/p/cosmos 170-190/65</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="2">quasar® lt med</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="3">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20	model:	quasar® lt med		U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 1600 VA / momentary 3500 VA			max weight on running surface: 300 kg / 660 lbs			<table><tr><td>(21)cos30003va19-0001 (11)151113</td><td> 2015-11-13</td><td></td></tr><tr><td> </td><td> h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td rowspan="2"></td></tr><tr><td>(01)4050588000516</td><td> cos30003va19-0001 Made in Germany</td></tr></table>	(21)cos30003va19-0001 (11)151113	 2015-11-13		 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588000516	 cos30003va19-0001 Made in Germany
product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20																									
model:	quasar® lt med																												
U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A																											
power input: long time 1600 VA / momentary 3500 VA																													
max weight on running surface: 300 kg / 660 lbs																													
(21)cos30003va19-0001 (11)151113	 2015-11-13																												
 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																												
(01)4050588000516	 cos30003va19-0001 Made in Germany																												
quasar® med cos30003va20	<table><tr><td>product family:</td><td colspan="2">treadmill h/p/cosmos 170-190/65</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="2">quasar® med</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="3">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20	model:	quasar® med		U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 1600 VA / momentary 3500 VA			max weight on running surface: 300 kg / 660 lbs			<table><tr><td>(21)cos30003va20-0001 (11)070928</td><td> 2007-09-28</td><td></td></tr><tr><td> </td><td> h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td rowspan="2"></td></tr><tr><td>(01)4050588000523</td><td> cos30003va20-0001 Made in Germany</td></tr></table>	(21)cos30003va20-0001 (11)070928	 2007-09-28		 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588000523	 cos30003va20-0001 Made in Germany
product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20																									
model:	quasar® med																												
U: 230 V ~	f: 50 Hz / 60 Hz	class: S, I, A																											
power input: long time 1600 VA / momentary 3500 VA																													
max weight on running surface: 300 kg / 660 lbs																													
(21)cos30003va20-0001 (11)070928	 2007-09-28																												
 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																												
(01)4050588000523	 cos30003va20-0001 Made in Germany																												
quasar® med 3p cos30003va26	<table><tr><td>product family:</td><td colspan="2">treadmill h/p/cosmos 170-190/65</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="2">quasar® med 3p</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="3">power input: long time 2700 VA / momentary 5700 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20	model:	quasar® med 3p		U: 400 V 3~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 2700 VA / momentary 5700 VA			max weight on running surface: 300 kg / 660 lbs			<table><tr><td>(21)cos30003va26-0001 (11)151111</td><td> 2015-11-11</td><td></td></tr><tr><td> </td><td> h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td rowspan="2"></td></tr><tr><td>(01)4050588001544</td><td> cos30003va26-0001 Made in Germany</td></tr></table>	(21)cos30003va26-0001 (11)151111	 2015-11-11		 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001544	 cos30003va26-0001 Made in Germany
product family:	treadmill h/p/cosmos 170-190/65		CE 0123 MD	IP20																									
model:	quasar® med 3p																												
U: 400 V 3~	f: 50 Hz / 60 Hz	class: S, I, A																											
power input: long time 2700 VA / momentary 5700 VA																													
max weight on running surface: 300 kg / 660 lbs																													
(21)cos30003va26-0001 (11)151111	 2015-11-11																												
 	 h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																												
(01)4050588001544	 cos30003va26-0001 Made in Germany																												

pulsar® lt cos30004va01	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>pulsar® lt</td><td>MD</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="2">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="2">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	pulsar® lt	MD	U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A			power input: long time 1600 VA / momentary 3500 VA		max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30004va01-0001 (11)111013</td><td></td><td>2011-10-13</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30004va01-0001</td></tr><tr><td>(01)4050588000554</td><td colspan="3">Made in Germany </td></tr></table>	(21)cos30004va01-0001 (11)111013		2011-10-13		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30004va01-0001	(01)4050588000554	Made in Germany 		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																														
model:	pulsar® lt	MD																															
U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A																																
power input: long time 1600 VA / momentary 3500 VA																																	
max weight on running surface: 300 kg / 660 lbs																																	
(21)cos30004va01-0001 (11)111013		2011-10-13																															
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
			cos30004va01-0001																														
(01)4050588000554	Made in Germany 																																
pulsar® lt 3p cos30004va02	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>pulsar® lt 3p</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz / 60 Hz class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="2">power input: long time 2700 VA / momentary 5700 VA</td></tr><tr><td colspan="2">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	pulsar® lt 3p	MD	U: 400 V 3~	f: 50 Hz / 60 Hz class: S, I, A			power input: long time 2700 VA / momentary 5700 VA		max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30004va02-0001 (11)151113</td><td></td><td>2015-11-13</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30004va02-0001</td></tr><tr><td>(01)4050588000561</td><td colspan="3">Made in Germany </td></tr></table>	(21)cos30004va02-0001 (11)151113		2015-11-13		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30004va02-0001	(01)4050588000561	Made in Germany 		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																														
model:	pulsar® lt 3p	MD																															
U: 400 V 3~	f: 50 Hz / 60 Hz class: S, I, A																																
power input: long time 2700 VA / momentary 5700 VA																																	
max weight on running surface: 300 kg / 660 lbs																																	
(21)cos30004va02-0001 (11)151113		2015-11-13																															
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
			cos30004va02-0001																														
(01)4050588000561	Made in Germany 																																
pulsar® cos30004va03	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>pulsar®</td><td>MD</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="2">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="2">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	pulsar®	MD	U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A			power input: long time 1600 VA / momentary 3500 VA		max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30004va03-0001 (11)071227</td><td></td><td>2007-12-27</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30004va03-0001</td></tr><tr><td>(01)4050588000578</td><td colspan="3">Made in Germany </td></tr></table>	(21)cos30004va03-0001 (11)071227		2007-12-27		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30004va03-0001	(01)4050588000578	Made in Germany 		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																														
model:	pulsar®	MD																															
U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A																																
power input: long time 1600 VA / momentary 3500 VA																																	
max weight on running surface: 300 kg / 660 lbs																																	
(21)cos30004va03-0001 (11)071227		2007-12-27																															
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
			cos30004va03-0001																														
(01)4050588000578	Made in Germany 																																
pulsar® 3p cos30004va04	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>pulsar® 3p</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz / 60 Hz class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="2">power input: long time 2700 VA / momentary 5700 VA</td></tr><tr><td colspan="2">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	pulsar® 3p	MD	U: 400 V 3~	f: 50 Hz / 60 Hz class: S, I, A			power input: long time 2700 VA / momentary 5700 VA		max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30004va04-0001 (11)071203</td><td></td><td>2007-12-03</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30004va04-0001</td></tr><tr><td>(01)4050588000585</td><td colspan="3">Made in Germany </td></tr></table>	(21)cos30004va04-0001 (11)071203		2007-12-03		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30004va04-0001	(01)4050588000585	Made in Germany 		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																														
model:	pulsar® 3p	MD																															
U: 400 V 3~	f: 50 Hz / 60 Hz class: S, I, A																																
power input: long time 2700 VA / momentary 5700 VA																																	
max weight on running surface: 300 kg / 660 lbs																																	
(21)cos30004va04-0001 (11)071203		2007-12-03																															
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
			cos30004va04-0001																														
(01)4050588000585	Made in Germany 																																
locomotion® 190/65 DE med cos30024va03	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>locomotion® 190/65 DE med</td><td>MD</td></tr><tr><td>U: 230 V ~</td><td>f: 50 Hz / 60 Hz class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="2">power input: long time 1600 VA / momentary 3500 VA</td></tr><tr><td colspan="2">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	locomotion® 190/65 DE med	MD	U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A			power input: long time 1600 VA / momentary 3500 VA		max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30024va03-0001 (11)151113</td><td></td><td>2015-11-13</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30024va03-0001</td></tr><tr><td>(01)4050588002107</td><td colspan="3">Made in Germany </td></tr></table>	(21)cos30024va03-0001 (11)151113		2015-11-13		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30024va03-0001	(01)4050588002107	Made in Germany 		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																														
model:	locomotion® 190/65 DE med	MD																															
U: 230 V ~	f: 50 Hz / 60 Hz class: S, I, A																																
power input: long time 1600 VA / momentary 3500 VA																																	
max weight on running surface: 300 kg / 660 lbs																																	
(21)cos30024va03-0001 (11)151113		2015-11-13																															
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
			cos30024va03-0001																														
(01)4050588002107	Made in Germany 																																

locomotion® 190/65-3p DE med cos30024va04	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 170-190/65</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>locomotion® 190/65-3p DE med</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz / 60 Hz</td><td>class: S, I, A</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="3">power input: long time 2700 VA / momentary 5700 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20	model:	locomotion® 190/65-3p DE med	MD	U: 400 V 3~	f: 50 Hz / 60 Hz	class: S, I, A			power input: long time 2700 VA / momentary 5700 VA			max weight on running surface: 300 kg / 660 lbs					<table><tr><td>(21)cos30024va04-0001 (11)111013</td><td></td><td>2011-10-13</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30024va04-0001 </td></tr><tr><td>(01)4050588002114</td><td colspan="3">Made in Germany</td></tr></table>	(21)cos30024va04-0001 (11)111013		2011-10-13		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30024va04-0001 	(01)4050588002114	Made in Germany		
product family:	treadmill h/p/cosmos 170-190/65	CE 0123	IP20																																	
model:	locomotion® 190/65-3p DE med	MD																																		
U: 400 V 3~	f: 50 Hz / 60 Hz	class: S, I, A																																		
power input: long time 2700 VA / momentary 5700 VA																																				
max weight on running surface: 300 kg / 660 lbs																																				
(21)cos30024va04-0001 (11)111013		2011-10-13																																		
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																		
			cos30024va04-0001 																																	
(01)4050588002114	Made in Germany																																			
venus® 200/75 cos30005-01va05	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>venus® 200/75</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A code: 3.1</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="3">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20	model:	venus® 200/75	MD	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1			power input: long time 7500 VA / momentary 15200 VA			max weight on running surface: 300 kg / 660 lbs					<table><tr><td>(21)cos30005-01va05-0001 (11)120124</td><td></td><td>2012-01-24</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30005-01va05-0001 </td></tr><tr><td>(01)4050588001612</td><td colspan="3">Made in Germany</td></tr></table>	(21)cos30005-01va05-0001 (11)120124		2012-01-24		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30005-01va05-0001 	(01)4050588001612	Made in Germany		
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20																																	
model:	venus® 200/75	MD																																		
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1																																		
power input: long time 7500 VA / momentary 15200 VA																																				
max weight on running surface: 300 kg / 660 lbs																																				
(21)cos30005-01va05-0001 (11)120124		2012-01-24																																		
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																		
			cos30005-01va05-0001 																																	
(01)4050588001612	Made in Germany																																			
venus® 200/75r cos30005-01va06	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>venus® 200/75r</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A code: 3.1</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="3">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20	model:	venus® 200/75r	MD	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1			power input: long time 7500 VA / momentary 15200 VA			max weight on running surface: 300 kg / 660 lbs					<table><tr><td>(21)cos30005-01va06-0001 (11)101220</td><td></td><td>2010-12-20</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30005-01va06-0001 </td></tr><tr><td>(01)4050588001629</td><td colspan="3">Made in Germany</td></tr></table>	(21)cos30005-01va06-0001 (11)101220		2010-12-20		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30005-01va06-0001 	(01)4050588001629	Made in Germany		
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20																																	
model:	venus® 200/75r	MD																																		
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1																																		
power input: long time 7500 VA / momentary 15200 VA																																				
max weight on running surface: 300 kg / 660 lbs																																				
(21)cos30005-01va06-0001 (11)101220		2010-12-20																																		
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																		
			cos30005-01va06-0001 																																	
(01)4050588001629	Made in Germany																																			
venus® 200/100 cos30006-01va05	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>venus® 200/100</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A code: 3.1</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="3">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20	model:	venus® 200/100	MD	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1			power input: long time 7500 VA / momentary 15200 VA			max weight on running surface: 300 kg / 660 lbs					<table><tr><td>(21)cos30006-01va05-0001 (11)180726</td><td></td><td>2018-07-26</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30006-01va05-0001 </td></tr><tr><td>(01)4050588001674</td><td colspan="3">Made in Germany</td></tr></table>	(21)cos30006-01va05-0001 (11)180726		2018-07-26		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30006-01va05-0001 	(01)4050588001674	Made in Germany		
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20																																	
model:	venus® 200/100	MD																																		
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1																																		
power input: long time 7500 VA / momentary 15200 VA																																				
max weight on running surface: 300 kg / 660 lbs																																				
(21)cos30006-01va05-0001 (11)180726		2018-07-26																																		
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																		
			cos30006-01va05-0001 																																	
(01)4050588001674	Made in Germany																																			
venus® 200/100r cos30006-01va06	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td>CE 0123</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>venus® 200/100r</td><td>MD</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A code: 3.1</td><td rowspan="2"></td><td rowspan="2"></td></tr><tr><td colspan="3">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="3">max weight on running surface: 300 kg / 660 lbs</td><td colspan="2"></td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20	model:	venus® 200/100r	MD	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1			power input: long time 7500 VA / momentary 15200 VA			max weight on running surface: 300 kg / 660 lbs					<table><tr><td>(21)cos30006-01va06-0001 (11)120807</td><td></td><td>2012-08-07</td><td></td></tr><tr><td rowspan="2"> </td><td rowspan="2"></td><td colspan="2">h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td></td><td>cos30006-01va06-0001 </td></tr><tr><td>(01)4050588001681</td><td colspan="3">Made in Germany</td></tr></table>	(21)cos30006-01va06-0001 (11)120807		2012-08-07		 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com			cos30006-01va06-0001 	(01)4050588001681	Made in Germany		
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123	IP20																																	
model:	venus® 200/100r	MD																																		
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A code: 3.1																																		
power input: long time 7500 VA / momentary 15200 VA																																				
max weight on running surface: 300 kg / 660 lbs																																				
(21)cos30006-01va06-0001 (11)120807		2012-08-07																																		
 		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																		
			cos30006-01va06-0001 																																	
(01)4050588001681	Made in Germany																																			

saturn® 250/75 cos30007-01va05	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="3">saturn® 250/75</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20	model:	saturn® 250/75			U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5			power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30007-01va05-0001 (11)111216</td><td>2011-12-16</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td></td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>(01)4050588001735</td><td>SN</td><td>cos30007-01va05-0001</td></tr><tr><td colspan="2"></td><td>Made in Germany</td></tr></table>	(21)cos30007-01va05-0001 (11)111216	2011-12-16	h/p/cosmos	UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	(01)4050588001735	SN	cos30007-01va05-0001			Made in Germany
product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20																																	
model:	saturn® 250/75																																					
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																																			
power input: long time 7500 VA / momentary 15200 VA																																						
max weight on running surface: 300 kg / 660 lbs																																						
(21)cos30007-01va05-0001 (11)111216	2011-12-16	h/p/cosmos																																				
UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																				
(01)4050588001735	SN	cos30007-01va05-0001																																				
		Made in Germany																																				
saturn® 250/75r cos30007-01va06	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="3">saturn® 250/75r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20	model:	saturn® 250/75r			U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5			power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30007-01va06-0001 (11)120521</td><td>2012-05-21</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td></td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>(01)4050588001742</td><td>SN</td><td>cos30007-01va06-0001</td></tr><tr><td colspan="2"></td><td>Made in Germany</td></tr></table>	(21)cos30007-01va06-0001 (11)120521	2012-05-21	h/p/cosmos	UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	(01)4050588001742	SN	cos30007-01va06-0001			Made in Germany
product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20																																	
model:	saturn® 250/75r																																					
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																																			
power input: long time 7500 VA / momentary 15200 VA																																						
max weight on running surface: 300 kg / 660 lbs																																						
(21)cos30007-01va06-0001 (11)120521	2012-05-21	h/p/cosmos																																				
UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																				
(01)4050588001742	SN	cos30007-01va06-0001																																				
		Made in Germany																																				
saturn® 250/100 cos30008-01va05	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="3">saturn® 250/100</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20	model:	saturn® 250/100			U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5			power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30008-01va05-0001 (11)171005</td><td>2017-10-05</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td></td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>(01)4050588001797</td><td>SN</td><td>cos30008-01va05-0001</td></tr><tr><td colspan="2"></td><td>Made in Germany</td></tr></table>	(21)cos30008-01va05-0001 (11)171005	2017-10-05	h/p/cosmos	UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	(01)4050588001797	SN	cos30008-01va05-0001			Made in Germany
product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20																																	
model:	saturn® 250/100																																					
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																																			
power input: long time 7500 VA / momentary 15200 VA																																						
max weight on running surface: 300 kg / 660 lbs																																						
(21)cos30008-01va05-0001 (11)171005	2017-10-05	h/p/cosmos																																				
UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																				
(01)4050588001797	SN	cos30008-01va05-0001																																				
		Made in Germany																																				
saturn® 250/100r cos30008-01va06	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="3">saturn® 250/100r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20	model:	saturn® 250/100r			U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5			power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30008-01va06-0001 (11)110804</td><td>2011-08-04</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td></td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>(01)4050588001803</td><td>SN</td><td>cos30008-01va06-0001</td></tr><tr><td colspan="2"></td><td>Made in Germany</td></tr></table>	(21)cos30008-01va06-0001 (11)110804	2011-08-04	h/p/cosmos	UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	(01)4050588001803	SN	cos30008-01va06-0001			Made in Germany
product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20																																	
model:	saturn® 250/100r																																					
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																																			
power input: long time 7500 VA / momentary 15200 VA																																						
max weight on running surface: 300 kg / 660 lbs																																						
(21)cos30008-01va06-0001 (11)110804	2011-08-04	h/p/cosmos																																				
UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																				
(01)4050588001803	SN	cos30008-01va06-0001																																				
		Made in Germany																																				
saturn® 250/125r cos30009-01va03	<table><tr><td>product family:</td><td colspan="3">treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td colspan="3">saturn® 250/125r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td><td rowspan="3"></td><td rowspan="3"></td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20	model:	saturn® 250/125r			U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5			power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30009-01va03-0001 (11)120417</td><td>2012-04-17</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td></td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td></tr><tr><td>(01)4050588001834</td><td>SN</td><td>cos30009-01va03-0001</td></tr><tr><td colspan="2"></td><td>Made in Germany</td></tr></table>	(21)cos30009-01va03-0001 (11)120417	2012-04-17	h/p/cosmos	UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	(01)4050588001834	SN	cos30009-01va03-0001			Made in Germany
product family:	treadmill h/p/cosmos 200-300/75-125			CE 0123 MD	IP20																																	
model:	saturn® 250/125r																																					
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																																			
power input: long time 7500 VA / momentary 15200 VA																																						
max weight on running surface: 300 kg / 660 lbs																																						
(21)cos30009-01va03-0001 (11)120417	2012-04-17	h/p/cosmos																																				
UDI		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																																				
(01)4050588001834	SN	cos30009-01va03-0001																																				
		Made in Germany																																				

saturn® 300/75 cos30010-01va05	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>saturn® 300/75</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20	model:	saturn® 300/75	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5	power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30010-01va05-0001 (11)201006</td><td>2020-10-06</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td></td></tr><tr><td>(01)4050588001889</td><td>SN cos30010-01va05-0001</td><td></td></tr><tr><td></td><td>Made in Germany</td><td></td></tr></table>	(21)cos30010-01va05-0001 (11)201006	2020-10-06	h/p/cosmos	UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001889	SN cos30010-01va05-0001			Made in Germany	
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20																													
model:	saturn® 300/75																															
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																													
power input: long time 7500 VA / momentary 15200 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30010-01va05-0001 (11)201006	2020-10-06	h/p/cosmos																														
UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
(01)4050588001889	SN cos30010-01va05-0001																															
	Made in Germany																															
saturn® 300/75r cos30010-01va06	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>saturn® 300/75r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20	model:	saturn® 300/75r	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5	power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30010-01va06-0001 (11)201006</td><td>2020-10-06</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td></td></tr><tr><td>(01)4050588001896</td><td>SN cos30010-01va06-0001</td><td></td></tr><tr><td></td><td>Made in Germany</td><td></td></tr></table>	(21)cos30010-01va06-0001 (11)201006	2020-10-06	h/p/cosmos	UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001896	SN cos30010-01va06-0001			Made in Germany	
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20																													
model:	saturn® 300/75r																															
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																													
power input: long time 7500 VA / momentary 15200 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30010-01va06-0001 (11)201006	2020-10-06	h/p/cosmos																														
UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
(01)4050588001896	SN cos30010-01va06-0001																															
	Made in Germany																															
saturn® 300/100 cos30011-01va05	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>saturn® 300/100</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20	model:	saturn® 300/100	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5	power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30011-01va05-0001 (11)201006</td><td>2020-10-06</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td></td></tr><tr><td>(01)4050588001940</td><td>SN cos30011-01va05-0001</td><td></td></tr><tr><td></td><td>Made in Germany</td><td></td></tr></table>	(21)cos30011-01va05-0001 (11)201006	2020-10-06	h/p/cosmos	UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001940	SN cos30011-01va05-0001			Made in Germany	
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20																													
model:	saturn® 300/100																															
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																													
power input: long time 7500 VA / momentary 15200 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30011-01va05-0001 (11)201006	2020-10-06	h/p/cosmos																														
UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
(01)4050588001940	SN cos30011-01va05-0001																															
	Made in Germany																															
saturn® 300/100r cos30011-01va06	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>saturn® 300/100r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20	model:	saturn® 300/100r	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5	power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30011-01va06-0001 (11)121016</td><td>2012-10-16</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td></td></tr><tr><td>(01)4050588001957</td><td>SN cos30011-01va06-0001</td><td></td></tr><tr><td></td><td>Made in Germany</td><td></td></tr></table>	(21)cos30011-01va06-0001 (11)121016	2012-10-16	h/p/cosmos	UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001957	SN cos30011-01va06-0001			Made in Germany	
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20																													
model:	saturn® 300/100r																															
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																													
power input: long time 7500 VA / momentary 15200 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30011-01va06-0001 (11)121016	2012-10-16	h/p/cosmos																														
UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
(01)4050588001957	SN cos30011-01va06-0001																															
	Made in Germany																															
saturn® 300/125r cos30012-01va03	<table><tr><td>product family:</td><td>treadmill h/p/cosmos 200-300/75-125</td><td rowspan="2">CE 0123 MD</td><td rowspan="2">IP20</td></tr><tr><td>model:</td><td>saturn® 300/125r</td></tr><tr><td>U: 400 V 3~</td><td>f: 50 Hz - 60 Hz</td><td>class: S, I, A</td><td>code: 4.5</td></tr><tr><td colspan="4">power input: long time 7500 VA / momentary 15200 VA</td></tr><tr><td colspan="4">max weight on running surface: 300 kg / 660 lbs</td></tr></table>	product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20	model:	saturn® 300/125r	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5	power input: long time 7500 VA / momentary 15200 VA				max weight on running surface: 300 kg / 660 lbs				<table><tr><td>(21)cos30012-01va03-0001 (11)110913</td><td>2011-09-13</td><td>h/p/cosmos</td></tr><tr><td>UDI</td><td>h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com</td><td></td></tr><tr><td>(01)4050588001988</td><td>SN cos30012-01va03-0001</td><td></td></tr><tr><td></td><td>Made in Germany</td><td></td></tr></table>	(21)cos30012-01va03-0001 (11)110913	2011-09-13	h/p/cosmos	UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com		(01)4050588001988	SN cos30012-01va03-0001			Made in Germany	
product family:	treadmill h/p/cosmos 200-300/75-125	CE 0123 MD	IP20																													
model:	saturn® 300/125r																															
U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5																													
power input: long time 7500 VA / momentary 15200 VA																																
max weight on running surface: 300 kg / 660 lbs																																
(21)cos30012-01va03-0001 (11)110913	2011-09-13	h/p/cosmos																														
UDI	h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com																															
(01)4050588001988	SN cos30012-01va03-0001																															
	Made in Germany																															

saturn® 450/300rs cos30013-02va02	product family: treadmill h/p/cosmos 450-300				CE 0123 	IP00
	model: saturn® 450/300rs					
	U: 400 V 3~	f: 50 Hz - 60 Hz	class: S, I, A	code: 4.5		
	current input: long time 21.0 A / momentary 44.0 A					
	max weight on running surface: 300 kg / 660 lbs					

(21)cos30013-02va02-0001		2022-03-14	h/p/cosmos
(11)220314		h/p/cosmos sports & medical gmbh 83365 Nussdorf-Traunstein / Germany service@hpcosmos.com	
			
(01)4050588002923		cos30013-02va02-0001	
		Made in Germany	