







	<h2 style="margin: 0;">order form</h2> <h3 style="margin: 0;">consumables & service</h3> <p style="margin: 0;">fax: +49 86 69 86 42 49</p> <p style="margin: 0;">email: service@hpcosmos.com</p>	
--	--	--

We hereby confirm our order for delivery on: _____ dispatch type: _____ urgent economic

pcs.	order number	description	picture	€ per unit	total €
	cos000200008001	Special oil (bottle 0.25 l) For h/p/cosmos running deck surface only for <u>running applications</u> (not for bicycle and wheel-chair use!) and <u>manual</u> lubrication with "hand-pump" syringe)		18.90 €	
	cos14007	Special oil (bottle 0.25 l) for automatic lubrication This oil type (special lubricant) is only for running deck surfaces for applications with wheels (cycles, wheelchair, ski, skating, etc.) where the running machine has an oil tank and <u>automatic</u> oil pump.		25.60 €	
	cos102446	Special oil for lubrication man.125ml for pluto treadmills		18.20 €	
	cos12181	Syringe 10ml For lubrication of running surface		1.10 €	
	cos10905	POLAR chest belt XS	 chest-belt rubber-band (only rubber band, without transmitter)	7.55 €	
	cos10906	POLAR chest belt S		7.55 €	
	cos10165	POLAR chest belt M		7.55 €	
	cos10907	POLAR chest belt L		7.55 €	
	cos10902	POLAR Transmitter Set T31 uncoded, with chest-belt rubber-band size M		41.70 €	
	cos15178	POLAR Transmitter Set T34 uncoded, <u>extended range</u> for running machines, with chest-belt rubber-band size M		55.00 €	
	cos14903-03-XXS	Chest belt system XXS size XXS (children) for h/p/cosmos safety arch		395.00 €	
	cos14903-03-XS	Chest belt system XS size XS (children) for h/p/cosmos safety arch		395.00 €	
	cos14903-03-S	Chest belt system S size S for h/p/cosmos safety arch		355.00 €	
	cos14903-03-M	Chest belt system M size M for h/p/cosmos safety arch		355.00 €	
	cos14903-03-L	Chest belt system L size L for h/p/cosmos safety arch		355.00 €	
	cos14903-03-XL	Chest belt system XL size XL for h/p/cosmos safety arch		395.00 €	
	cos13752	airwalk, cable assembly for h/p/cosmos airwalk 35/50/90	To be replaced every 6 months or earlier at the first sign of wear.	134.00 €	
	cos60000110029	1 time maintenance + safety check treadmill series mercury / mercury med	model:	on request. ask local dealer	
	cos60000110030	1 time maintenance + safety check treadmill series quasar/med ... pulsar 3p	serial-no.	on request	
	cos60000110032	1 time maintenance + safety check treadmill series venus/saturn up to 300/125	<small>Costs do not include any spare parts. Scope see general terms of business.</small>	on request	
	cos11964-V4.05.1-xxx	MCU4 EPROM Firmware V4.05.1 treadmill For treadmill with MCU4 control board with coscom v3		45.40 €	
	cos16600m5tm-v1.08.3	MCU5 Firmware v1.08.3 running machine For treadmill with MCU5 control board with coscom v3 and coscom v4		38.80 €	
Total price in € (only valid for Germany!), excl. packing and shipment costs, EXW, without VAT					€
Institution / Company & Delivery address:			Institution / Company & Billing address (if different):		
Place, Date:			Stamp and Signature of Customer		

Prices valid only for Germany and valid from January 1 – December 31, 2020. Prices in other countries can vary significantly. Transport, packing, VAT, import taxes, custom duties, L/C and bank fees, installation and instruction are not included. Full terms of trading available on request. E & OE. Subject to alterations without prior notice: http://www.h-p-cosmos.com/en/company/terms_of_business.htm

